

VistaLab  
Catalog #

Model

Color

Accuracy\*

Precision\*

1065	1mL pipette	Gray	±0.75%	0.50%
1066	2mL pipette	Red	±0.75%	0.50%
1067	3mL pipette	Green	±0.75%	0.50%
1068	4mL pipette	White	±0.75%	0.50%
1069	5mL pipette	Yellow	±0.75%	0.50%
1076	6mL pipette	Red	±0.60%	0.30%
1077	7mL pipette	Red	±0.60%	0.30%
1078	8mL pipette	Red	±0.60%	0.30%
1079	9mL pipette	Red	±0.60%	0.30%
1070	10mL pipette	Red	±0.60%	0.30%
1135	2,3,4,5mL pipette	Blue	±1.0%	0.50%
1135	2,3,5,10mL pipette	Red	±0.75%	0.50%
1061	1-5mL pipette	Blue	±1.0% at 5mL ±1.0% at 3mL ±1.0% at 1mL	±.50% at 5mL ±.50% at 3mL ±.60% at 1mL

